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Bruker AXS Microanalysis GmbH

QUANTAX

Energy-dispersive X-ray spectrometer
for electron microscopy

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



The screenshot displays the Bruker software interface for point analysis. The interface is divided into several sections:


- Top Bar:** Shows 'Point Mode' and various device parameters.
 - Devices:** Includes 'Spectrometer' (Input: 300 kcps, Range: 10 keV / 60 kcps, Dead time: 0 %, Real time: 9.7 kcps), 'Spectrometer' (Range: 20 keV / 60 kcps), and 'Imaging system' (Format: 800 x 600, HV: 15.0 kV, Magnification: 5766).
- Project:** Displays a list of projects, including 'Demo spectra'.
- Image scan:** Shows a grayscale SEM image of a sample. A green crosshair indicates the analysis point (labeled 3). The image includes technical details: 'SE MAG: 5766 x', 'HV: 15.0 kV', 'VWD: 20.0 mm', 'Continuously: 800 x 600', '28 x 21 μm'.
- Spectrum:** A graph showing 'cps/eV' vs energy. Peaks are labeled with elements: N, C, O, Sn, Cu, and Sn. The 'Acquire' button is highlighted (labeled 6).
- Measuring time:** A panel on the right with settings for 'Automatic' (Precise), 'Manual', 'Real time [s]', 'Live time [s]', 'Counts', 'Region start [keV]', and 'Region end [keV]'. It also includes 'Multiple spectrometer' and 'Cyclic acquisition' options.
- Automatic analysis:** A panel on the right with options for 'None', 'Continuous', and 'After acquisition' (selected).
- Coating correction:** A panel on the right with options for 'None', 'Coating correction Carbon only', and 'Coating calibration'.
- Spectrum name:** A panel on the right with fields for 'Acquisition', 'Spectrum number: 1', and 'Automatic numbering'.

Numbered callouts (1-7) indicate the steps for performing a point analysis:

1. Select 'Point' in the 'Project' section.
2. Press 'Preview' in the 'Image scan' section.
3. Adjust the analysis point position on the SEM image.
4. Adjust SEM beam parameters (e.g., Input, Range, Dead time) in the 'Devices' section.
5. Select 'Automatic' in the 'Measuring time' panel.
6. Press 'Acquire' in the 'Spectrum' section.
7. View the resulting spectrum graph.

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The screenshot displays the Bruker software interface for energy calibration. The interface is divided into several sections:

- System Navigation:** At the top, there are tabs for System, User interface, Spectrometer, Imaging, SEM, and Stage. The Spectrometer tab is selected.
- Energy calibration Spectrometer:** This section contains the main calibration settings.
 - Calibration settings:** Includes Pulse throughput (60000 cps, 90000 cps, 120000 cps, 275000 cps) and Energy range (10 keV, 20 keV, 40 keV, 80 keV). The 'Current' button is highlighted.
 - Calibration:** Includes Coarse, Mean, and Precise options. The 'Start' button is highlighted.
- Current values:** Displays Zeropeak channel (191.770), Channel width [eV] (2.501), FWHM at 0.0 keV (43.08), and FWHM at 8.0463 keV (144.45). A 'Test' button is present.
- Periodic Table:** A periodic table is shown with the element Copper (Cu) highlighted in red.
- Spectrum Plot:** A plot showing cps/eV vs keV. The x-axis ranges from 0 to 10 keV, and the y-axis ranges from 0 to 160 cps/eV. A prominent peak is labeled 'Cu' at 8.0463 keV.
- Confirm Dialog:** A dialog box asks 'Do you want to save new calibration data?' with 'Yes' and 'No' buttons.

Numbered callouts (1-6) indicate the steps for system calibration:

1. System → Spectrometer
2. Select parameters, press current to get current parameters
3. Select calibration element
4. Select calibration accuracy
5. Press start to start calibration
6. Press yes to confirm updating the database

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